Protect Your Patients

KYRA™ Armboards with optional KYRA CLEAN™ Armboard Pads

KYRA™ Gravity Lock Armboard

- Easy to use the weighted trigger arm board pivots to the desired position and locks in place.
- Easy on the hands ergonomic trigger allow the arm to be moved to any position on the table side-rail or to be removed.
- Measures 26" x 6" (66 x 15.25 cm).
- Patient weight capacity 500 lbs.

KYRA2000 KYRA Gravity Lock Armboard (US) KYRA2005 KYRA Gravity Lock Armboard (EU)









KYRA CLEAN™ Armboard Pad

- Seamless and free of VELCRO® brand fasteners* the KYRA CLEAN™ Armboard is easy to clean and designed to improve infection control.
- Easy to mount with no hook and loop fasteners.
- Soft, pressure relieving foam pad available.
- Available in 2" (5 cm) and 3" (7.5 cm) thickness.
- Fits armboards 26" x 6" (66 x 15.25 cm). (AMSCO style)

KYRA2020 KYRA CLEAN™ Armboard Pad 2" KYRA2025 KYRA CLEAN™ Armboard Pad 3"



KYRA COMFORT™ Armboard Pad

- Traditional pads with soft pressure relieving foam.
- Available in 2" (5 cm) and 3" (7.5 cm) thickness.
- Fits armboards 26" x 6" (66 x 15.25 cm).

KYRA2010 KYRA COMFORT™ Armboard Pad 2" KYRA2015 KYRA COMFORT™ Armboard Pad 3"





KYRA Clamps









KYRA SECURE-RELEASE™ Clamps

Mounts anywhere on the table rail in one easy motion. Requires no awkward hand motions.

- Simply push the clamp against the side rail to place it.
- Push the release button to release the clamps.
- Mounts any blade accessory.

KYRA4100 KYRA SECURE-RELEASE™ Clamp, US KYRA4200 KYRA SECURE-RELEASE™ Clamp, EU KYRA4300 KYRA SECURE-RELEASE™ Clamp, UK/ES KYRA4400 KYRA SECURE-RELEASE™ Clamp, Denyer







Autoclave Safe KYRA RAIL™ Clamp

- Mounts anywhere on the table rail for fast and easy mounting.
- Rounded ergonomic shape is easy on the hands.
- Mounts any blade accessory.
- Autoclavable

KYRA4600 KYRA RAIL™ Clamp, US KYRA4610 KYRA RAIL™ Clamp, EU KYRA4620 KYRA RAIL™ Clamp, UK/ES KYRA4630 KYRA RAIL™ Clamp, Denyer









innovative solutions for today's modern OR™ 888-611-KYRA (5972) | kyramedical.com